

FOR RESEARCH USE ONLY!

DiaEasy™ Dialyzer (3 ml) Floating racks

(Catalog # 1016-5, -10; Store at RT [15-25°C])

8/13

I. CATALOG #	1016-5	5 Floating racks
	1016-10	10 Floating racks

II. DESCRIPTION:

Biovision's DiaEasy™ dialyzer (3 ml) floating racks are designed to support our DiaEasy™ Dialyzer tubes of 3 ml volume. Each of these racks can hold upto 4 tubes.

III. STORAGE AND HANDLING:

Store the racks at room temperature.

IV. COMPANION PRODUCTS:

- DiaEasy™ Dialyzer (10, 15, 20 ml) Floating racks (Cat #. 1000-5, -10)
- DiaEasy™ Dialyzer (3 ml) Supporting trays (Cat #. 1017-5, -10)
- DiaEasy™ Dialyzer (800 µl) Floating racks (Cat #. 1018-5, -10)
- DiaEasy™ Dialyzer (800 µl) Supporting trays (Cat #. 1019-5, -10)
- DiaEasy™ Dialyzer (250 µl) Floating racks (Cat #. 1020-5, -10)
- DiaEasy™ Dialyzer (250 µl) Supporting trays (Cat #. 1028-5, -10)




Figure 1: DiaEasy™ Dialyzer (3 ml) Floating rack.

V. Related Products

DiaEasy™ Dialyzer (250 µl) MWCO 6-8 kDa (K1020-10, -25, -100)	DiaEasy™ Dialyzer (250 µl) MWCO 12-14 kDa (K1021-10, -25, -100)
DiaEasy™ Dialyzer (800 µl) MWCO 3.5 kDa (K1018-10, -25, -100)	DiaEasy™ Dialyzer (800 µl) MWCO 6-8 kDa (K1019-10, -25, -100)
DiaEasy™ Dialyzer (3 ml) MWCO 50 kDa (K1016-10, -25, -100)	DiaEasy™ Dialyzer (800 µl) MWCO 1 kDa (K1017-10, -25, -100)
DiaEasy™ Dialyzer (3 ml) MWCO 3.5 kDa (K1012-10, -25, -100)	DiaEasy™ Dialyzer (3 ml) MWCO 6-8 kDa (K1013-10, -25, -100)
DiaEasy™ Dialyzer (3 ml) MWCO 12-14 kDa (K1014-10, -25, -100)	DiaEasy™ Dialyzer (3 ml) MWCO 25 kDa (K1015-10, -25, -100)
DiaEasy™ Dialyzer (20 ml) MWCO 3.5 kDa (K1000-10, -25, -100)	DiaEasy™ Dialyzer (15 ml) MWCO 3.5 kDa (K1001-10, -25, -100)
DiaEasy™ Dialyzer (10 ml) MWCO 3.5 kDa (K1002-10, -25, -100)	DiaEasy™ Dialyzer (20 ml) MWCO 6-8 kDa (K1003-10, -25, -100)
DiaEasy™ Dialyzer (15 ml) MWCO 6-8 kDa (K1004-10, -25, -100)	DiaEasy™ Dialyzer (10 ml) MWCO 6-8 kDa (K1005-10, -25, -100)
DiaEasy™ Dialyzer (20 ml) MWCO 12-14 kDa (K1006-10, -25, -100)	DiaEasy™ Dialyzer (15 ml) MWCO 12-14 kDa (K1007-10, -25, -100)
DiaEasy™ Dialyzer (10 ml) MWCO 12-14 kDa (K1008-10, -25, -100)	DiaEasy™ Dialyzer (20 ml) MWCO 1 kDa (K1009-10, -25, -100)
DiaEasy™ Dialyzer (15 ml) MWCO 1 kDa (K1010-10, -25, -100)	

FOR RESEARCH USE ONLY! Not to be used on humans.