

Apolipoprotein A4, Human Plasma

CATALOG #: 7282-50 50 µg

ALTERNATE NAMES: Apolipoprotein A-IV, apo AIV, Apo-AIV, ApoA-IV

SOURCE: Human plasma. Prepared from plasma shown to be non-reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.

PURITY: ≥ 85% by SDS - PAGE

MOL. WEIGHT: 46 kDa


FORM: Frozen

FORMULATION: In 1 mM NH₄HCO₃, pH 8.2

STORAGE CONDITIONS: Aliquot and store at -70 °C or lower. Avoid repeated freezing and thawing cycles.

DESCRIPTION:

Apolipoprotein AIV (Apo AIV) is a glycoprotein synthesized by the human intestine. Synthesis appears to be stimulated by the absorption of high density lipoproteins that are rich in triglycerides rather than by the uptake or re-esterification of fatty acids. The formation of chylomicrons acts as a signal for the induction of Apo AIV synthesis a signal enhanced by a factor from the ileum, probably PYY. Apo AIV alters the activity of the key enzymes (LPL and LCAT) of lipoprotein metabolism and cholesterol efflux from extra hepatic tissues. Also, since stimulation of intestinal synthesis and secretion by lipid absorption are rapid, Apo AIV likely plays a role in the short-term regulation of food intake. Other evidence suggests involvement in the long-term regulation of both food intake and body weight. Chronic ingestion of high-fat diet blunts the intestinal Apo AIV response to lipid feeding and, consequently, may explain why the chronic ingestion of a high-fat diet predisposes both animals and humans to obesity. It is required for efficient activation of lipoprotein lipase by ApoC-II and is a potent activator of LCAT. Apo AIV is an amphipathic protein that can emulsify lipids and has been linked to protective roles against cardiovascular disease and obesity.



Apolipoprotein AIV, Human Plasma

4-12% Bis-Tris-HCl NuPAGE gel
 1. MW Standard
 2. Apo AIV - 5 µg (reduced / heated)
 3. Apo AIV - 10 µg (reduced / heated)
 4. Apo AIV - 5 µg (non-reduced / no heat)
 5. Apo AIV - 10 µg (non-reduced / no heat)

RELATED PRODUCTS:

- ApoA-1, human recombinant (Cat # 4693-20, -100, -1000)
- ApoA-1, rat recombinant (Cat # 4695-20, -100, -1000)
- Apo-D, human recombinant (Cat # 4326-10, -1000)
- APO-J / Clusterin (Cat. No. 4325-10, -500)
- Apo-SAA, human recombinant (Cat # 4324-50, -500, -1000)
- ApoE2, human recombinant (Cat # 4760-50, -100, -500, -1000, -5000)
- ApoE3, human recombinant (Cat # 4696-50, -100, -500, -5000)
- ApoE4, human recombinant (Cat # 4699-50, -100, -500, -1000, -5000)
- Apolipoprotein CI, Human Plasma (Cat # 4703-100)
- Apolipoprotein CII, Human Plasma (Cat # 4704-50)
- Apolipoprotein CIII, Human Plasma (Cat # 4706-100)

FOR RESEARCH USE ONLY! Not to be used in humans.