

Caspase-Family Inhibitor, Z-VAD-FMK

CATALOG #:	1010-20C
LOT #:	_____
STORAGE CONDITIONS:	Store at -20°C.
SHELF LIFE:	1 year under proper storage conditions
MOLECULAR WEIGHT:	468.0
SEQUENCE:	Z-Val-Ala-Asp-FMK
PURITY:	One spot by TLC analysis
PACKAGE CONTENTS:	20 µl (10 mM) of Z-VAD-FMK in DMSO

DESCRIPTION:

A synthetic peptide that irreversibly inhibits activity of caspase family proteases and blocks apoptosis. The inhibitor is designed as a methyl ester to facilitate cell permeability. The product can be used in both in vitro and in vivo studies.

RECOMMENDED USAGE:

We recommend using 1000X dilutions for inhibiting caspase activity in Jurkat cell culture (e.g., Add 1 µl to 1 ml of culture medium). The optimal doses may vary for different cells and culture conditions.

NOTE: For research purpose only! Not to be used in humans.




Fig. 1. Apoptosis was induced in Jurkat cells with camptothecin (2 µM) for 6 hours in the presence of various concentrations of Caspase Family Inhibitor, Z-VAD-FMK. Caspase activity was measured using Caspase-3 Fluorometric Assay Kit (K105-25) according to the kit instructions.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction Molecules