

Caspase-8/FLICE Inhibitor, Z-IETD-FMK

CATALOG #:	1064-20C
LOT #:	_____
STORAGE CONDITIONS:	Store at -20°C.
SHELF LIFE:	1 year under proper storage conditions
MOLECULAR WEIGHT:	654.0
SEQUENCE:	Z-Ile-Glu-Thr-Asp-FMK (FMK, fluoromethyl ketone)
PURITY:	>99% by HPLC analysis.
PACKAGE CONTENTS:	20 Z-IETD-FMK (10 mM) in DMSO

DESCRIPTION: A synthetic peptide that irreversibly inhibits FLICE and related protease/caspase activity and blocks apoptosis. IETD is the specific recognition sequence found in CPP32 proenzyme. The inhibitor is designed as a methyl ester to facilitate cell permeability.

RECOMMENDED USAGE:

We recommend using 5000X dilutions for inhibiting caspase-8 activity in Jurkat cell culture (e.g., Add 1 | to 5 ml of culture medium). The optimal doses may vary for different cells and culture conditions.

REFERENCES:

1. Thornberry, N.A., *et al.* (1997) *J. Biol.Chem.* **272**:17907-17911.
2. Nicholson, D.W., *et al.* (1995) *Nature* **376**:37-40.

RELATED PRODUCTS:

1. KinaseSTAR™ JNK & Akt Activity Assay Kit
2. Histone Deacetylase Fluorometric & Colorimetric Assay Kits
3. Calpain Activity Assay Kit
4. Apoptotic Cell Isolation Kit
5. MitoCapture™ Mitochondrial Apoptosis Detection Kit
6. Mitochondria/Cytosol Fractionation Kit
7. Nuclear/Cytosol Fractionation Kit
8. Cytochrome c Releasing Apoptosis Assay Kit
9. ApoSENSOR™ ATP Determination Kit & ADP/ATP Ratio Assay Kit
10. ApoGSH™ Glutathione Detection Kit
11. Nitric Oxide Fluorometric & Colorimetric Detection Kits
12. Annexin V-FITC, -Cy3, -Cy5, -EGFP, -PE, -Biotin Kits and Bulk Reagents
13. CaspSCREEN™ Flow Cytometric Caspase Assay Kit
14. CaspGLOW™ Active Caspase Staining Kits
15. CaspSELECT™ Caspase-3 & -7 Immunoassay Kits
16. Caspase-1,-2,-3,-5,-6,-8,-9,-10 Activity Assay Kits
17. Recombinant Human, Mouse, Rat Active Caspases & Procaspsases
18. Ready-to-use Caspase Inhibitors & Sets
19. Ready-to-use Caspase Substrates and Sets
20. Quick Apoptotic DNA Ladder Detection Kit
21. TUNEL-based DNA Fragmentation Assay Kits
22. Ready-to-use Apoptosis Inducers and Set
23. Quick Cell Proliferation Assay Kit
24. LDH-Cytotoxicity Assay Kit
25. Live-Dead Cell Staining Kit
26. Anti-Caspase-1,-2,-3,-4,-5,-6,-7,-8,-9,-10,-11,-12,-13,-14 & Active Caspases
27. Anti-Apaf-1, Bcl-2, Bcl-x, Bax, Bad, Bid, Bim/Bod, Mcl-1, Cytochrome c
28. Anti-DFF40, DFF45, FADD, p53, IRAK, Wee1, CD40, DcR2, cFLIP, CTGF
29. Many Other Apoptosis & Cell Signaling Related Antibodies
30. Anti-GFP, His, GST, and Other Tag Antibodies
31. Growth Factors, Cytokines & Chemokines (many)
32. Link-FAST™ 5 Minutes DNA Ligation Kit
32. Gel-FAST™ 20 Minutes Gel Staining/Destaining Kit
33. Protease Inhibitor Cocktail