

<b>CATALOG #:</b>	1108-20C
<b>LOT #:</b>	_____
<b>STORAGE CONDITIONS:</b>	Store at -20°C.
<b>SHELF LIFE:</b>	6 months under proper storage conditions
<b>MOLECULAR WEIGHT:</b>	652.0
<b>SEQUENCE:</b>	Z-Leu-Glu-Val-Asp-FMK (FMK, fluoromethyl ketone)
<b>PURITY:</b>	>95% by HPLC analysis.
<b>PACKAGE CONTENTS:</b>	20 µl (10 mM) Z-LEVD-FMK in DMSO

### DESCRIPTION:

A synthetic peptide inhibitor that inhibits caspase-4 and related caspase activity. FMK inhibitors are irreversible and non-toxic, and can be used in both *in vivo* and *in vitro* studies.

### RECOMMENDED USAGE:

We recommend using 1000-5000X dilutions for inhibiting caspase-4 and related caspase activity (e.g., Add 1 µl to 1-5 ml of culture medium). However, the optimal doses may vary for different cells and culture conditions. The inhibitor is designed as a methyl ester to facilitate cell permeability.

### REFERENCES:

1. Thornberry, N.A., *et al.* (1997) *J. Biol. Chem.* **272**:17907-17911.

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit

