


PRODUCT: Phorbol 12-myristate 13-acetate (PMA)

CATALOG #: 1544-5

AMOUNT: 5 mg

MOLECULAR FORMULA: $C_{36}H_{56}O_8$

MOLECULAR WEIGHT: 616.8

CAS NUMBER: 16561-29-8

PURITY: $\geq 99\%$

SOLUBILITY: Soluble in DMSO or ethanol

STORAGE CONDITIONS: Store at -20°C . Protect from light.

DESCRIPTION: This is the most commonly used phorbol ester. It binds to and activates protein kinase C, causes an extremely wide range of effects in cells and tissues and is a very potent mouse skin tumor promoter. Inhibits apoptosis induced by the Fas antigen, but induced apoptosis in HL-60 promyelocytic leukemia cells.

CAUTION: Highly irritant tumor promoter

FOR RESEARCH USE ONLY! Not to be used in humans.**RELATED PRODUCTS:**

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

