

PRODUCT: Ionomycin, free acid

CATALOG #:	1565-5
AMOUNT:	5 mg
LOT #:	_____
MOLECULAR FORMULA:	C ₄₁ H ₇₂ O ₉
MOLECULAR WEIGHT:	709.02
CAS:	[56092-81-01]
PURITY:	≥97%
APPEARANCE:	Off-white to light yellow Powder
SOLUBILITY:	Soluble in DMSO
STORAGE CONDITIONS:	-20°C in the dark. Do not expose to sunlight or fluorescent lights for extended periods.
BACKGROUND:	Antibiotic that acts as a potent and selective calcium ionophore; more effective than A23187. Ionomycin binds Ca ²⁺ in the 7.0-9.5 pH range with the resulting complex exhibiting intense UV absorption. Ionomycin induces apoptosis in immature B cell lines and in cultured embryonic rat cortical neurons. More mature B cell lines do not undergo apoptosis, but are arrested in the G1 phase of the cell cycle.
CAUTION:	Tumor promoter. Do not ingest, swallow, or inhale. Do not get in eyes or on skin. Care should be taken to prevent contact through all routes of exposure.

For Research Use Only! Not to be used in humans.

FEATURED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers & Inhibitors
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit

Growth Factors and Cytokines (many)

Polyclonal and Monoclonal Antibodies