

PRODUCT: 10Kd Spin Column

CATALOG #:	1997-25
QUANTITY:	25
MEMBRANE:	10,000 MWCO Polyethersulfone
CONCENTRATOR CUP:	Polycarbonate
FILTRATE VESSEL:	Polypropylene
STORAGE:	Room Temperature

DESCRIPTION:

The 10Kd Spin Columns are disposable ultrafiltration devices for the concentration of biological samples. Maximum initial sample volume 500 μ l. They can be used in either a swing bucket or fixed angle rotors accepting 2.0 ml centrifuge tube at maximum speed 15000 g. The dead stop volume is ~5-15 μ l.

The vertical membrane design and thin channel filtration chamber minimizes membrane fouling and provides high speed concentration, even laden solutions. It is specifically designed with low internal surface and membrane area in order to achieve superior recovery from dilute solution.



Cat# 1997-25 (Contains 25 10Kd Spin Columns)

FOR RESEARCH USE ONLY! Not to be used in human.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies