

## 3X SDS-PAGE Loading Buffer

**CATALOG #:** 2108-10

**AMOUNT:** 5 x 2 ml

**LOT #:** \_\_\_\_\_

**STORAGE CONDITIONS:** Short term: Store at 4°C.  
Long term: Store at -20°C.

**SHELF LIFE:** Stable for 1 year at -20°C

**BUFFER CONTENTS:**

Component	Concentration
Tris-HCl (pH 6.8)	150 mM
DTT	300 mM
SDS	6%
Bromophenol blue	0.3%
Glycerol	30%

**For Research Use Only! Not to be used in human.**

**RELATED PRODUCTS:**

- Apoptosis Detection Kits & Reagents
  - Annexin V Kits & Bulk Reagents
  - Caspase Assay Kits & Reagents
  - Mitochondrial Apoptosis Kits & Reagents
  - Nuclear Apoptosis Kits & Reagents
  - Apoptosis Inducers and Set
  - Apoptosis siRNA Vectors
- Cell Fractionation System
  - Mitochondria/Cytosol Fractionation Kit
  - Nuclear/Cytosol Fractionation Kit
  - Membrane Protein Extraction Kit
  - Cytosol/Particulate Rapid Separation Kit
  - Mammalian Cell Extraction Kit
  - FractionPREP Fractionation System
- Cell Proliferation & Senescence
  - Quick Cell Proliferation Assay Kit
  - Senescence Detection Kit
  - High Throughput Apoptosis/Cell Viability Assay Kits
  - LDH-Cytotoxicity Assay Kit
  - Bioluminescence Cytotoxicity Assay Kit
  - Live/Dead Cell Staining Kit
- Cell Damage & Repair
  - HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
  - HAT Colorimetric Assay Kit & Reagents
  - DNA Damage Quantification Kit
  - Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits
- Signal Transduction
  - cAMP & cGMP Assay Kits
  - Akt & JNK Activity Assay Kits
  - Beta-Secretase Activity Assay Kit
- Adipocyte & Lipid Transfer
  - Recombinant Adiponectin, Survivin, & Leptin
  - CETP Activity Assay & Drug Discovery Kits
  - PLTP Activity Assay & Drug Discovery Kits
  - Total Cholesterol Quantification Kit
- Molecular Biology & Reporter Assays
  - siRNA Vectors
  - Cloning Insert Quick Screening Kit
  - Mitochondrial & Genomic DNA Isolation Kits
  - 5 Minutes DNA Ligation Kit
  - 20 Minutes Gel Staining/Destaining Kit
  - $\beta$ -Galactosidase Staining Kit & Luciferase Reporter Assay Kit
- Growth Factors and Cytokines
- Monoclonal and Polyclonal Antibodies