

# Bad Polyclonal Antibody

**CATALOG #:** 3030-100

**AMOUNT:** 100 µg

**LOT #:** \_\_\_\_\_

**HOST:** Rabbit

**IMMUNOGEN:** Synthetic peptide surrounding amino acid 115 of mouse Bad

**FORMULATION:**

100 µg (0.5 mg/ml) antigen affinity purified rabbit anti-Bad polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA, 0.02% thimerosal.

**SPECIES REACTIVITY:** human, mouse, monkey, rat

**STORAGE CONDITIONS:**

Store at -20°C. For long term storage, aliquot and refreeze at -70°C. Avoid repeated freeze/thaw cycles.

**DESCRIPTION:**

The Bcl-2 (B-cell leukemia/lymphoma 2) family proteins function to either inhibit (Bcl-2 and Bcl-x<sub>L</sub>) or promote (Bax and Bad) apoptosis. Similar to other family member, Bad also exhibits sequence homology to Bcl-2. Bad can selectively dimerize with Bcl-x<sub>L</sub> and Bcl-2, but not with Bax, Bcl-x<sub>s</sub>, Mcl-1, or itself. In mammalian cells, Bad binds more strongly to Bcl-x<sub>L</sub> than Bcl-2. This may explain why Bad reverses the death repressor activity of Bcl-x<sub>L</sub>, but not that of Bcl-2

**SPECIFICITY:**

The antibody detects the 22 kDa Bad protein and it does not cross-react with related family members.

**APPLICATION AND USAGE:**

The antibody can be used for Western blot analysis (1-2 µg/ml), immunoprecipitation (4-8 µg/ml), and immunohistochemistry (20-40 µg/ml). However, the optimal conditions should be determined individually.

**For Research Use Only! Not to be used in human subjects.**

**RELATED PRODUCTS:**

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies