

JNK Polyclonal Antibody

CATALOG #:	3701-100
AMOUNT:	100 µg
LOT #:	_____
HOST:	Rabbit
IMMUNOGEN:	Synthetic peptide surrounding amino acid 412 of human JNK-1 (Internal ID# BV-520)
SPECIES REACTIVITY:	human, mouse, rat

FORMULATION:

100 µg (0.5 mg/ml) affinity purified rabbit anti-JNK polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 30% glycerol, 0.5% BSA, 0.01% thimerosal.

STORAGE CONDITIONS:


Store at -20°C. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

BACKGROUND:

The JNK pathway is activated by largely distinct stimuli including inflammatory cytokines (such as TNF- α and IL-1), UV light, inhibitors of protein synthesis and osmotic stress. Activated MEKK1 phosphorylates SEK1 (also known as MKK4), which in turn activates SAPK (also known as JNK). JNK binds tightly to the N-terminal region of c-Jun and ATF-2, and phosphorylates c-Jun at Ser63 and Ser73 and ATF-2 at Thr69 and Thr71.

APPLICATION AND USAGE:

The antibody can be used for Western blotting (0.5-4 µg/ml). Reacts with human, and in a lesser extent with mouse and rat. Per researchers feed back it can also be used in immunoprecipitation (10-20 µg/ml), and immunohistochemistry (10-20 µg/ml). However, the optimal conditions should be determined individually. Blocking peptide (Cat.# 3701BP-50) is available separately.



Cat #: 3701-100
Western blot analysis of JNK expression in Jurkat cell lysates (Lane 1, 2).

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies