

# Caper Polyclonal Antibody

<b>CATALOG #:</b>	3733-100
<b>AMOUNT:</b>	100 $\mu$ L
<b>LOT #:</b>	_____
<b>HOST:</b>	Rabbit
<b>IMMUNOGEN:</b>	Synthetic peptide of human Caper
<b>SPECIES REACTIVITY:</b>	human, mouse

## FORMULATION:

100  $\mu$ L affinity purified rabbit polyclonal antibody in phosphate-buffered saline (PBS) containing 30% glycerol, 0.5% BSA and 0.01% thimerosal.

## STORAGE CONDITIONS:

Store at  $-20^{\circ}\text{C}$ . For long term storage, aliquot and freeze at  $-70^{\circ}\text{C}$ . Avoid repeated freeze/thaw cycles.

## ALTERNATIVE NAMES:

RNPC2, RNA-binding region containing protein 2, Caper, Coactivator of activating protein-1 and estrogen receptors, HCC1, Hepatocellular carcinoma protein 1, Splicing factor CC1.3, Splicing factor CC1.4

## SPECIFICITY:

The antibody recognizes the Caper of human and mouse origins. Reactivity to other species has not been tested.

## APPLICATION:

The antibody can be used in Western blotting (1:500 – 1:2500) and Immunoprecipitation. HeLa cell lysate can be used as a positive control. However, the optimal concentrations should be determined individually.

**WARRANTY:** Warranty information can be found in [www.biovision.com](http://www.biovision.com).

**USAGE:** For Research Use Only. Not to be used in humans.

## RELATED PRODUCTS:

### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- $\beta$ -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

### Growth Factors and Cytokines

### Monoclonal and Polyclonal Antibodies