

CRSP1(TRAP220) Polyclonal Antibody

CATALOG #:	3735-100
AMOUNT:	100 µL
LOT #:	_____
HOST:	Rabbit
IMMUNOGEN:	Synthetic peptide between amino acid 1500 and 1600 of human CRSP1.

FORMULATION:

100 µl purified rabbit polyclonal antibody in phosphate-buffered saline (PBS) containing 30% glycerol, 1% BSA and 0.02% thimerosal.

SPECIES REACTIVITY: human, mouse

STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

ALTERNATIVE NAMES:

CRsp1, CRsp-1, Cofactor required for Sp1 Transcriptional Activation 1, TRAP220, Thyroid Hormone Receptor-Associated Protein Complex Component, 220 kD, PPARBP, PPAR Binding Protein, DRIP205, Vitamin D Receptor Interacting Protein Complex Component, 205 kD, DRIP230, vitamin D Receptor-Interacting Protein, 230 kD, TRIP2, Thyroid Receptor Interacting Protein 2, MED1.PBP, PPARGBP, RB18A.

SPECIFICITY:

The antibody recognizes the CRSP1 of human and mouse origins. Reactivity to other species has not been tested.

APPLICATION:

The antibody can be used in Western blotting (1:500 – 1:2500 µg/ml), Immunoprecipitation. However, the optimal concentrations should be determined individually.

USAGE: For Research Use Only. Not to be used in humans.

WARRANTY: Warranty information can be found in www.biovision.com.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies