

# MAD1 Polyclonal Antibody

**CATALOG #:** 3765-100

**AMOUNT:** 100  $\mu$ L

**LOT #:** \_\_\_\_\_

**HOST:** Rabbit

**IMMUNOGEN:** Synthetic peptide of human MAD1

**FORMULATION:**

100  $\mu$ l affinity purified rabbit polyclonal antibody in phosphate-buffered saline (PBS) containing 30% glycerol, 1% BSA and 0.02% thimerosal.

**SPECIES REACTIVITY:** human, mouse

**STORAGE CONDITIONS:**

Store at  $-20^{\circ}\text{C}$ . For long term storage, aliquot and freeze at  $-70^{\circ}\text{C}$ . Avoid repeated freeze/thaw cycles.

**ALTERNATIVE NAMES:**

MAD1, MAD-1, MADL1, HsMAS1, Mitotic arrest deficient-like 1, PIG9,p53, inducible protein 9, TP53I9,tumor protein p53 inducible protein 9, TXBP181, Tax bindingprotein-181

**SPECIFICITY:**

The antibody recognizes the MAD1 of human and mouse origins. Reactivity to other species has not been tested.

**APPLICATION:**

The antibody can be used in Western blotting (1:500 – 1:2000) and Immunoprecipitation. HeLa cell lysate can be used as a positive control. However, the optimal concentrations should be determined individually.

**WARRANTY:** Warranty information can be found in [www.biovision.com](http://www.biovision.com).

**USAGE:** For Research Use Only. Not to be used in humans.

**RELATED PRODUCTS:**

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- $\beta$ -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies