

## Recombinant Human HGF

<b>CATALOG #:</b>	4509-10	10 µg
	4509-1000	1 mg
<b>LOT #:</b>	_____	
<b>SOURCE:</b>	<i>Baculovirus infected High-5 cells</i>	
<b>PURITY:</b>	>98% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of HGF.	
<b>MOL. WEIGHT:</b>	78 kDa	
<b>FORM:</b>	Lyophilized from 0.1 M L-Arginine HCl + 10 mM Tris, pH 8.0	

### RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.5-1 mg/ml. This solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or – 20°C for future use.

### STORAGE CONDITIONS:

The lyophilized HGF is best-stored desiccated below 0°C. Reconstitute HGF should be stored in working aliquots at -20°C.

### DESCRIPTION:

Human Hepatocyte Growth Factor (HGF) is also known as Scatter Factor and Hepatopoeitin A. HGF is a potent mitogen for mature parenchymal hepatocyte cells and acts as a growth factor for a broad spectrum of tissues and cell types. Recombinant human HGF is an 80 kDa disulfide-linked heterodimeric protein consisting of the alpha chain (463 amino acids) and the beta chain (234 amino acids).

### BIOLOGICAL ACTIVITY:

The ED<sub>50</sub> determined by the scattering activity in the MDCK cell assay is typically at 0.5-1.0 ng/ml.

**FOR RESEARCH USE ONLY! Not to be used in humans.**

### FEATURED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Monoclonal and Polyclonal Antibodies