

Recombinant Human BMP-13

CATALOG #:

4639-10	10 µg
4639-50	50 µg
4639-1000	1 mg

LOT #: _____

SOURCE: *E. coli*

PURITY: >95% by SDS-PAGE
Endotoxin level is <0.1 ng per µg of BMP-13.

FORM: Lyophilized with no additives

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute to 0.1-1.0 mg/ml in water containing BSA (50 µg BSA per 1 µg of protein). Further dilutions should be done in water and stored at 4°C for 1 week or -20°C for future use. Note: Due to solubility reasons, the protein should be kept at low pH.

STORAGE CONDITIONS:

The lyophilized protein is best-stored desiccated below 0°C. Reconstituted BMP-13 should be stored in working aliquots at -20°C.

DESCRIPTION:

BMPs (bone morphogenetic proteins) belong to the TGF-β superfamily of structurally related signaling proteins. These ligands mediate numerous physiological processes. BMPs initiate, promote and maintain cartilage and bone development, growth, remodeling and repair, in both prenatal development and postnatal growth of eye, heart, kidney, skin, and other tissues. Expression of BMP-13 has been found in the hypertrophic chondrocytes of ossifying long bone centers. BioVision's recombinant human BMP-13/CDMP-2 is a 27 kDa homodimeric disulfide-linked protein consisting of two 120 amino acids polypeptide chains.

BIOLOGICAL ACTIVITY:

The ED₅₀ as determined by its ability to induce alkaline phosphatase production by ATDC-5 chondrogenic cells is in the range of 2-3 µg/ml.

USAGE: For In vitro Research Use Only. Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies