

Recombinant Hemagglutinin

CATALOG #:	4844-10
AMOUNT:	10 µg
LOT #:	30644
SOURCE:	<i>Insect cells</i>
PURITY:	>90% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of Hemagglutinin.
MOL. WEIGHT:	57.5 kDa
FORM:	Sterile liquid solution (10 µg/6 µl) containing 10 mM Sodium Phosphate, pH 7.2, 150 mM NaCl.
STORAGE CONDITIONS:	Centrifuge the vial prior to opening. This rHA solution is stable for 2 weeks at 15°C, but should be stored at 4°C. For long-term storage, it is recommended to add a carrier protein (0.1% HSA or BSA). Freeze/thaw should be avoided.
DESCRIPTION:	Hemagglutinin is glycosylated with N-linked sugars, produced using baculovirus vector in insect cells and purified using proprietary chromatographic techniques.

For laboratory research use only! Not to be used in humans.

FEATURED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies