

Anti-Human MIP-3 α Polyclonal Antibody

CATALOG #: 5491-100

AMOUNT: 100 μ g

LOT #: _____

HOST: Rabbit

IMMUNOGEN: *E. coli*-expressed human MIP-3 α

FORMULATION:

100 μ g (0.5 mg/ml) antigen affinity purified rabbit anti-human MIP-3 α polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, and 0.02% thimerosal.

SPECIES REACTIVITY: human

STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

BACKGROUND DESCRIPTION:

Human macrophage inflammatory protein-3 alpha (MIP-3 alpha, Exodus, LARC) is a member of the beta-chemokine (C-C) family of cytokines. Human MIP-3 alpha is a 8.0 kDa protein containing 70 amino acid residues.

APPLICATIONS:

The antibody can be used for Western blot analysis (0.5-1 μ g/ml) and ELISA (0.5-4 μ g/ml). However, the optimal conditions should be determined individually.

USAGE: For Research Use Only. Not to be used in human or diagnostic procedures.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies