


# ELISA Kits for Cancer Research

Enzyme-linked Immunosorbent Assay (ELISA) is a powerful analytics tool to quantify proteins and small molecules in serum, plasma, tissue homogenates and other biological fluids. BioVision offers fully validated, accurate, and ready-to-use ELISA Kits for convenient and sensitive quantitation of various cancer and tumor related targets.

## Key Features:

- **Accurate:** reproducible results with low intra & inter assay variability
- **High Sensitivity and Specificity**
- **Versatile:** broad range of sample types such as serum, plasma, tissue, cell culture supernatant, urine and many other biological fluids.



**Figure:** Standard curve for human PD-L1 ELISA Kits (Left) and human Interleukin-10 ELISA kits (Right).

Product Name	Cat. No.	Specificity	Detection Type
PD-1 ELISA Kit	K4153-100	Human	Sandwich
PD-L1 ELISA Kit	K4155-100	Human	Sandwich
AKT1/PKB ELISA Kit	K4166-100	Human	Sandwich
Angiogenin ELISA Kit	K4802-100	Human	Sandwich
Angiopoietin-1 ELISA Kit	K7115-100	Human	Sandwich
Angiopoietin-2 ELISA Kit	K7116-100	Human	Sandwich
Angiostatin ELISA Kit	K4755-100	Human	Sandwich
ANGPTL3 Serum ELISA Kit	K4914-100	Human	Sandwich
ANGPTL3 Serum ELISA Kit	K4915-100	Mouse, Rat	Sandwich
ANGPTL6 Serum ELISA Kit	K4916-100	Human	Sandwich
BDNF ELISA Kit	K4788-100	Human	Sandwich
Beta-NGF ELISA Kit	K4787-100	Human	Sandwich
BLC/BCA-1/CXCL13 ELISA Kit	K4766-100	Human	Sandwich

Product Name	Cat. No.	Specificity	Detection Type
BMP-4 ELISA Kit	K4745-100	Human	Sandwich
CA19-9 ELISA Kit	K7427-100	Human	Sandwich
Cancer Antigen 125 ELISA Kit	K4803-100	Human	Sandwich
Cancer Antigen 15-3 (CA15-3) ELISA Kit	K4804-100	Human	Sandwich
Carcinoembryonic Antigen (CEA) ELISA Kit	K4805-100	Human	Sandwich
Caspase-9 ELISA Kit	K4165-100	Human	Sandwich
Cathepsin B ELISA Kit	K7408-100	Human	Sandwich
Cathepsin D ELISA Kit	K7409-100	Human	Sandwich
Chorionic Gonadotropin (hCG) ELISA Kit	K7424-100	Human	Sandwich
Cortisol (human/mouse/rat) ELISA Kit	K7430-100	Human, Mouse, Rat	Sandwich
CTLA4 /CD152 ELISA Kit	K4157-100	Human	Sandwich
EGFR/ErbB ELISA Kit	K4168-100	Human	Sandwich
GM-CSF ELISA Kit	K4772-100	Human	Sandwich
Growth Hormone (human) ELISA Kit	K7412-100	Human	Sandwich
HGF ELISA Kit	K4781-100	Human	Sandwich
IL-10 ELISA Kit	K4167-100	Human	Sandwich
IL-1 $\alpha$ ELISA Kit (For Lysates)	K4792-100	Human	Sandwich
IL-1 $\beta$ ELISA Kit	K4794-100	Human	Sandwich
IL-1 $\beta$ ELISA Kit	K4795-100	Mouse	Sandwich
IL-4 ELISA Kit	K4164-100	Human	Sandwich
IL-5 ELISA Kit	K4156-100	Human	Sandwich
IL-7 ELISA Kit	K4154-100	Human	Sandwich
Leptin ELISA Kit	K4777-100	Human	Sandwich
p53 ELISA Kit	K4829-100	Human	Sandwich
Progranulin ELISA Kit	K4738-100	Human	Sandwich
Prostate Specific Antigen ELISA Kit	K7431-100	Human	Sandwich
TGF- $\alpha$ ELISA Kit	K4339-100	Human	Sandwich
TGF- $\beta$ 1 ELISA Kit	K4342-100	Human	Sandwich
TNF- $\alpha$ ELISA Kit	K4779-100	Human	Sandwich
TNF- $\alpha$ ELISA Kit	K1051-100	Mouse	Sandwich
TNF- $\beta$ ELISA Kit	K4780-100	Human	Sandwich
VEGF ELISA Kit	K5363-100	Human	Sandwich
3-Nitrotyrosine (3-NT) ELISA Kit	K4158-100	Universal	Competitive
8-hydroxy-2'-deoxyguanosine (8-OHdG) ELISA Kit	K4160-100	Universal	Competitive