

BioVision's Exosome Products

Exosome Products


Studying Exosome has never been So user friendly !

- ELISA Kits
- FACS Kits
- Isolation Tools: Immunobeads, Immunobeads
- Precipitation Reagents
- Size Exclusion Columns
- Lyophilized Exosomes
- DNA Extraction Kits
- RNA Extraction Kits
- Antibodies
- Exosome Labeling
- Fluorescent Exosome Standards

Exosomes are small endosome derived lipid nanoparticles (50-120 nm in diameter), actively secreted by exocytosis in most living cells.

Biovision is excited to offer a wide selection of products dedicated to Exosome Research. This includes high quality reagents and kits that are quick, affordable and easy to use. Our new line of exosome products will facilitate your understanding of exosomes using a wide variety of samples.

The full range of our exosome products are described in the catalogue and can be purchased online at <http://www.biovision.com/products/exosome-research.html>



Product Highlights:

- Easy, quick, efficient, cheap and convenient exosome isolation tools (ExoPure™) including Isolation Kits (filtration based); Isolation columns (size exclusion, chromatography based); Reagents (precipitation based); and Immunobeads or Immunoplates (immunocapture based)
- Ready to use exosome quantification ELISA assay kits (ExoQuant™)
- Ready to use exosome FACS kits (ExoFACS™)
- Efficient exosome capture and high quality RNA extraction kits (ExoRNA™)
- Efficient exosome capture and DNA purification kits (ExoDNA™)
- Validated multi-application exosomal monoclonal and polyclonal antibodies (ExoAb™)
- Highly purified, lyophilized exosome standards (ExoStd™) as calibration standards
- Fluorescent labeled exosome standards for EV tracking studies

Product Name	Cat. No.	Size
ExoAb™ Exosome Antibodies	A1500-A1514	50 µg
*ExoDNA™ Exosome DNA extraction kits	K1230-K1231	Varies
ExoDNA™ DNA Extraction kits	K1242-K1244	Varies

155 S. Milipitas Blvd, Milipitas, CA 95035

T: 408-493-1800 F: 408-493-1801

Toll Free: 800-891-9699 (US Only)



Product Name	Cat. No.	Size
ExoQuant™ Exosome ELISA Assay Kits	K1201-K1210	100 Assays
*ExoFACS™ Exosome FACS Assay Kits	K1232-K1236	20 Reactions
ExoPure™ Exosome Isolation Columns	K1246-K1247	Varies
ExoPure™ Exosome Isolation Immunobeads	M1030-M1039	Varies
ExoPure™ Exosome Isolation Immunoplates	M1004-M1027	100 Assays
ExoPure™ Exosome Isolation Kits	K1237- K1241	Varies
*ExoPure™ Exosome Isolation Reagent	M1001-M1003	Varies
ExoRNA™ Exosome RNA Extraction kits	K1220-K1222	Varies
ExoStd™ Fluorescent Exosome Standard	M1073-M1079; M1083-M1091	Varies
ExoStd™ Lyophilized Exosome Standards	M1040-M1071	Varies

**Not available for sale in USA*

We are constantly adding new products in this field to provide a complete product portfolio for your exosome research. Please visit our website, <http://www.biovision.com/products/exosome-research.html> for the most up-to-date details on all the products that we have in our catalog.

For ordering please email us at: orders@biovision.com