


Fluorescent Proteins

Fluorescent Proteins (EGFP/RFP/YFP/BFP/CFP) were expressed and purified from transformed E. coli using a method that ensures high purity and maximal fluorescence. Fluorescent proteins are suitable as control reagents for expression studies. Applications include: Use as standards for SDS-PAGE, Western blot analysis, calibration of fluorometers and flow cytometers, fluorescence microscopy, and microinjection into cells and tissues, etc. The recombinant fluorescent proteins can also be conjugated to other proteins. Antibodies against these fluorescent proteins, mammalian expression vectors, as well as a fluorescence protein quantitation kit are also available.



EGFP - Green Fluorescent Protein

Product Name	Cat. No.	Size
GFP Polyclonal Antibody	3999	100 µg
GFP Quantitation Kit	K815	100 assays
pVision-GFP-C Vector	9998	20 µg
pVision-GFP-N Vector	9999	20 µg
Recombinant EGFP (GREEN)	4999	100 µg, 1 mg, 5 mg

RFP - Red Fluorescent Protein

Product Name	Cat. No.	Size
dsRED Polyclonal Antibody	3993	100 µg
mCherry Polyclonal Antibody	5993	100 µg
pVision-RFP-C Vector	9996	20 µg
pVision-RFP-N Vector	9997	20 µg
Recombinant mCherry RFP	4993	100 µg, 1 mg, 5 mg
Recombinant RFP (dsRED)	4997	100 µg, 1 mg, 5 mg

YFP - Yellow Fluorescent Protein

Product Name	Cat. No.	Size
Recombinant YFP (YELLOW)	4998	100 µg, 1 mg, 5 mg
YFP Polyclonal Antibody	3991	100 µg

155 S. Milipitas Blvd, Milipitas, CA 95035

T: 408-493-1800 F: 408-493-1801

Toll Free: 800-891-9699 (US Only)



CFP - Cyan Fluorescent Protein


Product Name	Cat. No.	Size
CFP Polyclonal Antibody	5986	100 µg
Recombinant CFP (CYAN)	4996	100 µg, 1 mg, 5 mg

BFP - Blue Fluorescent Protein

Product Name	Cat. No.	Size
BFP Polyclonal Antibody	5994	100 µg
Recombinant BFP (BLUE)	4994	100 µg, 1 mg, 5 mg

Related Sets

Product Name	Cat. No.	Size
Fluorescent Protein Antibodies Set	K817	6 x 100 µg
Fluorescent Protein Set	K816	6 x 100 µg



Please visit www.biovision.com for a comprehensive overview on our products.