

Lipid Metabolism - Products at a glance

BioVision offers an array of tools for lipid metabolism research – assay kits, proteins and antibodies to name a few. Below are some selected products. The comprehensive list can be found at:

<http://www.biovision.com/diabetes-obesity-metabolic-syndrome-1607/lipids-metabolism-1647/>

Assay Kits

Product Name	Cat #	Size	Product Name	Cat #	Size
3T3-L1 Differentiation Kit	K579-100	For 100 ml differentiation medium	PicoProbe™ β -Hydroxybutyrate Assay Kit (F)	K651-100	100 assays
Acetoacetate Assay Kit (C)	K650-100	100 assays	PicoProbe™ Lipolysis (3T3-L1) Assay Kit (F)	K578-100	100 assays
Adipogenesis Assay Kit (C/F)	K610-100	100 assays	PicoProbe™ Triglyceride (TG) Assay Kit (F)	K614-100	100 assays
CETP Activity Assay Kit (F)	K601-100	100 assays	PLTP Activity Assay Kit (F)	K604-100	100 assays
CETP Activity Assay Kit (F) II	K595-100	100 assays	PLTP Activity Assay Kit (F) II	K593-100	100 assays
CETP Inhibitor Drug Screening Kit (F)	K602-100	100 assays	PLTP Inhibitor Drug Screening Kit (F)	K605-100	100 assays
CETP Inhibitor Screening Kit (F) II	K594-100	100 assays	Preadipocyte Isolation Kit	K583-5	For 5 g Tissue
Cholesterol Detection Kit (cell-based)	K587-100	100 assays	Sphingomyelin Quantification Assay Kit (C)	K600-100	100 assays
Cholesterol Efflux Assay Kit (F) (cell-based)	K582-100	100 assays	Sphingomyelinase Activity Assay Kit (C)	K599-100	100 assays
Choline/Acetylcholine Quantification Kit (C/F)	K615-100	100 assays	Sphingomyelinase Activity Assay Kit (F)	K574-100	100 assays
Free Fatty Acid Quantification Kit (C/F)	K612-100	100 assays	Total Cholesterol and Cholesteryl Ester Assay Kit (C/F)	K603-100	100 assays
Free Glycerol Assay Kit (C) II	K634-100	100 assays	Triglyceride Quantification Kit (C/F)	K622-100	100 assays
Free Glycerol Assay Kit (C/F)	K630-100	100 assays			
Glucose Assay Kit (C/F)	K606-100	100 assays			
Glycerol-3-Phosphate Dehydrogenase Activity Assay Kit (C)	K640-100	100 assays			
HDL and LDL/VLDL Quantification Kit (C/F)	K613-100	100 assays			
HDL Uptake Assay Kit (F)	K586-100	100 assays			
HMG-CoA Reductase Activity/Inhibitor Screening Kit (C)	K588-100	100 assays			
β -Hydroxybutyrate (β -HB) Assay Kit (C)	K632-100	100 assays			
LDL Uptake Assay Kit (F)	K585-100	100 assays			
Lipase Activity Assay Kit (C)	K722-100	100 assays			
Lipolysis (3T3-L1) Assay Kit (C)	K577-100	100 assays			
Lipolysis (Adipocyte) Assay Kit (C/F)	K581-5	For 5 g Tissue			
Lipoprotein Lipase Activity Assay Kit (F)	K721-100	100 assays			
PAF Acetylhydrolase (PAF-AH) Inhibitor Screening Kit (C)	K766-100	100 assays			
PAF Acetylhydrolase Activity Assay Kit (C)	K765-100	100 assays			
Phosphatidylcholine Assay Kit (C/F)	K576-100	100 assays			
PhosphoSeek™ Sphingosine Kinase-1 Assay Kit (F)	K708-400	100 assays			
PicoProbe™ Free Glycerol Assay Kit (F)	K643-100	100 assays			

Antibodies

Product Name	Cat #	Size
ACADS Antibody	3539-100	100 μ g
ACADSB Antibody	3150-100	100 μ g
ACADSB Antibody	3157-100	100 μ g
ACADVL Antibody	3240-100	100 μ g
ACAT1 Antibody	3245-100	100 μ g
ACAT2 Antibody	3254-100	100 μ g
Adiponectin Antibody	5901-250	250 μ g
Adiponectin Antibody	5901-50	50 μ g
Adiponectin Antibody	5902-250	250 μ g
Adiponectin Antibody	5902-50	50 μ g
Adiponectin Antibody	5903-250	250 μ g
Adiponectin Antibody	5903-50	50 μ g
Adiponectin Antibody (Clone MADI-1147)	5900-50	50 μ g
ApoA1 Antibody	5693-100	100 μ g
ApoA1 Antibody	5695-100	100 μ g
ApoE3 Antibody	5696-100	100 μ g
ApoE4 Antibody	5699-100	100 μ g
ApoLipoprotein B Antibody	3661-100	100 μ g

Antibodies cont...

Product Name	Cat #	Size
CETP Antibody	3413-100	100 µg
CETP Polyclonal Antibody	6606-200	200 µg
FABP1 Antibody	5501-200	200 µg
FABP3 Antibody	5503-100	100 µg
FABP4 Antibody	3717-100	100 µl
FABP4 Antibody	5504-100	100 µg
FABP5 Antibody	5505-100	100 µg
FABP7 Antibody	5507-100	100 µg
FATP4 Antibody	3268-100	100 µg
LDLR Antibody	3839-100	100 µg
LRP-5/6 Antibody	3714-100	100 µg
LRP-5/6 Antibody (Clone 1A12)	3801-100	100 µg
LXR Antibody	3641-100	100 µg

Proteins

Product Name	Cat #	Size
ACAT-2, hu rec	4910-20, 1000	20 µg, 1 mg
ApoA-1, hu rec	4693-20, 100, 1000	20 µg, 100 µg, 1 mg
ApoA-1, rt rec	4695-20, 100, 1000	20 µg, 100 µg, 1 mg
ApoE2, hu rec	4760-50, 500, 1000, 5000	50 µg, 500 µg, 1 mg, 5 mg
ApoE3, hu rec	4696-50, 500, 1000, 5000	50 µg, 500 µg, 1 mg, 5 mg
ApoE4, hu rec	4699-50, 500, 1000, 5000	50 µg, 500 µg, 1 mg, 5 mg
Apolipoprotein All, human plasma	7281-500	500 µg
Apolipoprotein AIV, human plasma	7282-50	50 µg
Apolipoprotein B, human plasma	7283-500	500 µg
Apolipoprotein E, human plasma	7284-50	50 µg
Chylomicrons, human plasma	7285-1000	1 mg
CRP, hu rec	4863-10, 50, 1000	10 µg, 50 µg, 1 mg
CRP, hu rec	4864-250, 1000	250 µg, 1 mg
FABP1, hu rec	4501-25, 100, 1000	25 µg, 100 µg, 1 mg
FABP3, hu rec	4503-25, 100, 1000	25 µg, 100 µg, 1 mg
FABP4, hu rec	4504-25, 100, 1000	25 µg, 100 µg, 1 mg
FABP5, hu rec	4505-25, 100, 1000	25 µg, 100 µg, 1 mg
FABP7, hu rec	4507-25, 100, 1000	25 µg, 100 µg, 1 mg
Human CellExp™ FABP-3, hu rec	7234-10	10 µg
Human CellExp™ FABP2 /I-FABP, hu rec	7474-10, 50	10 µg, 50 µg
Carbonic Anhydrase 3, hu rec	7833-10, 50	10 µg, 50 µg
FNDC4, human/ms rec	7573-10	10 µg

Proteins cont...

Product Name	Cat #	Size
FNDC5, human/ms rec	7574-10	10 µg
FTO, hu rec	7557-10	10 µg
FTO, ms rec	7558-10	10 µg
Human CellExp™ Carbonic Anhydrase 10/CA10, hu rec	7485-10	10 µg
Human CellExp™ Carbonic Anhydrase 2/CA2, hu rec	7479-10, 50	10 µg, 50 µg
Human CellExp™ Carbonic Anhydrase 4/CA4, hu rec	7484-10	10 µg
Human CellExp™ Carbonic Anhydrase 9/CA9, hu rec	7478-10	10 µg
Human CellExp™ FOLR2, hu rec	7471-10, 50	10 µg, 50 µg
Human CellExp™ LDLR, hu rec	7472-10, 50	10 µg, 50 µg
Human CellExp™ Lipocalin-2 / LCN2, hu rec	7460-10, 50	10 µg, 50 µg
Human CellExp™ Lipocalin-2 / LCN2, ms rec	7568-10, 50	10 µg, 50 µg
Human CellExp™ Lipocalin-2 / LCN2, rt rec	7569-10, 50	10 µg, 50 µg
Human CellExp™ Meteorin-like, ms rec	7577-10	10 µg
Human CellExp™ Netrin-1, hu rec	7575-10	10 µg
Human CellExp™ RBP4, hu rec	7563-10, 50	10 µg, 50 µg
Human CellExp™ RBP4, ms rec	7564-10, 50	10 µg, 50 µg
Human CellExp™ RBP4, rt rec	7566-10, 50	10 µg, 50 µg
Human CellExp™ UNC5B, hu rec	7576-10	10 µg
Human CellExp™ VLDLR, hu rec	7464-10, 50	10 µg, 50 µg
RBP4, ms rec	7565-10, 50	10 µg, 50 µg
RBP4, rt rec	7567-10, 50	10 µg, 50 µg

Special discount:
Get 15% off ANY product
www.biovision.com
Code: VA15
Ends on 8/10/15

