

Molecular Biology Research Products

BioVision is delighted to offer a broad range of molecular biology products including DNA and RNA Purification Kits; Virus Purification Kits; DNA Clean Up Kits; PCR, QPCR, RT PCR, QRT PCR Master Mixes, Kits, Enzymes and Reagents; Molecular Biology Proteins and Enzymes; siRNA Expression Vectors; Luciferase Reporter Assay Kits and Substrates; DNA Ladders and Dyes and Numerous Molecular Biology Biochemicals including Antibiotics, Proteinase K, Transfection Reagents and many more. Over the years, our scientists have engineered the highest quality enzymes and most advanced and sophisticated formulations for PCR related applications. All of the BioVision's products have been well validated by researchers around the world, with an increasing number of citations available in scientific publications.

Product Highlights:

- Fast, efficient, safe and convenient DNA purification kits including mini, maxi, mega, midi and endotoxin free kits.
- Fast, convenient RNA purification kits from various samples including cells, tissues, bacteria and plants.
- Ready-to-use, fast, efficient, high recovery virus purification kits for Adenovirus, Adeno-associated Virus, Lentivirus, Retrovirus, and HCV.
- Ready-to-use, convenient and efficient DNA Clean Up Kits for recovering or concentrating a broad range of DNA fragments from agarose gel, PCR, or other enzymatic reactions.
- Ready-to-use, highly sensitive, efficient and cost-effective PCR, RT PCR, QPCR master mixes, kits, enzymes and reagents.
- Improved and efficient versions of proteins and enzymes such as CRISPR/Cas9 and many more.
- Molecular biology antibodies & blocking Peptides.
- siRNA expression vectors targeting important genes.
- DNA Ladders and Dyes.
- Molecular Biology Biochemicals including Antibiotics, Proteinase K, Luciferase substrates.
- High performance, animal-free, broad spectrum DNA and siRNA transfection reagents.
- Spin columns, test strips and many more products...

Some of the Key Molecular Biology Products:

Product Description	Cat. No.
DNA Purification Kits (mini, maxi, mega, midi and endotoxin free kits)	K1312-K1325; K1326-K1336 and More
RNA Purification Kits	K1340-K1358
Virus Concentration and Purification Kits (Adenovirus, Adeno-associated Virus, Lentivirus, Retrovirus, HCV)	K904; K1300-K1310
DNA Clean Up, DNA Quantification, DNA Ligation Kits and More	K1362-K1371; K900; M1217 and Many More
PCR, QPCR, RT PCR, QRT PCR Master Mixes, Kits, Enzymes and Reagents	M1105-M1190 and More
Molecular Biology Proteins/Enzymes; Antibodies and Blocking Peptides	7680, 7847 and Many More
Gene Expression Tools including Reporter Assay Kits, siRNA Expression Vectors, Activity Assays and More	9509 and Many More
DNA Ladders and Dyes	M1105-M1158; M1191-M1199
Molecular Biology Biochemicals (Antibiotics, Proteinase K, Transfection Reagents, Luciferase Substrates...)	1041; 9211 and Many More

<http://www.biovision.com/products/molecular-biology-tools.html>

We are constantly adding new products in this field to provide a complete product portfolio for your research. Please visit our website, www.biovision.com for the most up-to-date details on all the products that we have in our catalog.

For ordering please email us at: orders@biovision.com



Life Science SOURCE™
Innovation. Flexibility. Affordability

155 S Milpitas Blvd.
Milpitas, CA 95035
Tel: 408-493-1800
Fax: 408-493-1801