

Active Human Cathepsin H

CATALOG NO:	1023-5
AMOUNT:	5 µg
BACKGROUND:	A more basic protein than Cathepsin B or L. Functions as both an aminopeptidase and endopeptidase. Cathepsin H, like cathepsin L and B, is involved in the catabolism of proteins in the lysosomal system.
SOURCE:	Purified from human liver
FORM:	Frozen liquid (15 µl) in 50 mM sodium acetate, pH 5.5, with 1 mM EDTA and 100 mM NaCl
PURITY:	≥95% by SDS-PAGE
STORAGE CONDITIONS:	Store at -80°C, stable for at least 1 year as supplied. Avoid freeze/thaw cycles.
APPLICATIONS:	Human Cathepsin H is useful in inhibitor screening, kinetic studies, cleavage of substrates and other <i>in vitro</i> applications.
SPECIFIC ACTIVITY:	1.96 units/mg protein. One unit is defined as the amount of enzyme that hydrolyzes one µmol Arg-AFC substrate per minute at 37°C, pH 6.8, with 750 mM potassium phosphate, 1 mM EDTA and 3 mM cysteine.

Prepared from tissue shown to be non reactive for HbsAG, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.

For Research Use Only! Not to be used in human subjects.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit

Polyclonal and Monoclonal Antibodies

Growth Factors and Cytokines