

DiaEasy™ Dialyzer 96-Well Format Floating Rack

(Catalog # 1024-96-1; Store at RT)

04/21

I. **CATALOG #** 1024-96-1 **One 96-well format floating rack**

II. **DESCRIPTION:**

Biovision's DiaEasy™ dialyzer 96-Well Format Floating Rack is designed to support our DiaEasy™ Dialyzer tubes of 250 µl volume in a convenient 96-well format to accommodate multi-channel pipetting.

III. **STORAGE AND HANDLING:**

Store the rack at room temperature.

IV. **RELATED PRODUCTS:**

- DiaEasy™ Dialyzer (250 µl) Floating Racks (Cat. No. 1020)
- DiaEasy™ Dialyzer (250 µl) MWCO 6-8 kDa (Cat. No. K1020)
- DiaEasy™ Dialyzer (250 µl) MWCO 12-14 kDa (Cat. No. K1021)
- DiaEasy™ Dialyzer (250 µl) MWCO 25 kDa (Cat. No. K1022)
- DiaEasy™ Dialyzer (250 µl) Supporting Trays (Cat. No. 1028)




Figure 1: DiaEasy™ Dialyzer (250 µl) 96-Well Format Floating Rack, loaded with up to 96 DiaEasy™ Dialyzer tubes of 250 µl volume in the traditional 8 x 12 format.




Figure 2: DiaEasy™ Dialyzer (250 µl) 96-Well Format Floating Rack holding DiaEasy™ Dialyzer tubes of 250 µl volume during dialysis in the desired buffer.

FOR RESEARCH USE ONLY! Not to be used on humans.