

## Active Recombinant Human Caspase-4

**CATALOG NO:** 1084-25 25 units  
1084-100 100 units

### UNIT DEFINATION:

One unit of the recombinant caspase-4 is the enzyme activity that cleaves 1 nmol of the caspase substrate WEHD-pNA (pNA: p-nitroaniline) per hour at 37°C in a reaction solution containing 50 mM Hepes, pH 7.2, 50 mM NaCl, 0.1% Chaps, 10 mM EDTA, 5% Glycerol, and 10 mM DTT.

**FORM:** Lyophilized powder

**RECONSTITUTION:** Reconstitute to 1 unit per µl in distilled H<sub>2</sub>O

**PURITY:** ≥95% by SDS-PAGE

### STORAGE CONDITIONS:

The lyophilized caspase-4 is stable for 1 year at -70°C. Following reconstitution in water, the enzyme should be aliquoted and immediately stored at -70°C. Avoid multiple freeze/thaw cycles as activity might decrease.

### DESCRIPTION:

Caspase-4 (also know as Ich-4, ICE<sub>relII</sub>, TX) is a member of the caspase-family of cysteine proteases. Similar to other caspases, caspase-4 exists in cells as an inactive proenzyme which is matured by proteolysis to yield large and small subunits. The active caspase-4 is a heterotetramer consisting of two large and two small subunits. To date the mechanism of caspase-4 activation and roles of caspase-4 in apoptosis are poorly understood.

The recombinant active human caspase-4 was expressed in *E. coli*. The active caspase-4 is routinely tested at BioVision for its ability to enzymatically cleave these two substrates Ac-WEHD-pNA (Cat. #1102-200) or Ac-WEHD-AFC (1101-200).

### APPLICATIONS AND USAGE:

Active caspase-4 is useful in studying enzyme regulation, determining target substrates, screening caspase inhibitors, or as a positive control in caspase activity assays. We recommend using 1 unit/assay for analyzing caspase activity.

For a complete caspase-4 assay protocol, please refer to BioVision's Caspase-4 Fluorometric or Colorimetric Assay Kits (Cat.#: K126 and K127).

**FOR RESEARCH USE ONLY! Not to be used in humans.**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Quality Antibodies for Apoptosis and Signal Transduction Molecules