

## Active Recombinant Human Caspase-10/a

<b>CATALOG NO:</b>	1090A-25	25 units
	1090A-100	100 units

### UNIT DEFINITION:

One unit of the recombinant caspase-10/a is the enzyme activity that cleaves 1 nmol of the caspase substrate IETD-pNA (pNA: p-nitroanaline) per hour at 37° C in a reaction solution containing 50 mM Hepes, pH 7.2, 50 mM NaCl, 0.1% Chaps, 10 mM EDTA, 5% Glycerol, and 10 mM DTT.

**FORM:** 30% (w/v) Ammonium sulfate/water

**SPECIFIC ACTIVITY:** 5000 units/mg

**PURITY:** ≥90% by SDS-PAGE

### STORAGE CONDITIONS:

The caspase-10 is stable for 1 year at -70° C. Following reconstitution in PBS+15% glycerol, the enzyme should be aliquoted and immediately stored at -70° C. Avoid multiple freeze/thaw cycles as activity might decrease.

### DESCRIPTION:

Caspase-10/a (also known as Mch4) is a member of the caspase-family of cysteine proteases. Similar to other caspases, caspase-10 also exists in cells as an inactive proenzyme. During apoptosis procaspase-10 is processed at aspartate residues to form active caspase-10 which consists of two large (18 kDa) and two small (11 kDa) subunits in a heterotetramer form.

The recombinant active human caspase-10 was expressed in *E. coli*. The active caspase-10 is routinely tested at BioVision for its ability to enzymatically cleave these two substrates AEVD-AFC (Cat. No.1113-200) or AEVD-pNA (Cat. No. 1114-200).

### APPLICATIONS AND USAGE:

Active caspase10/a is useful in studying enzyme regulation, determining target substrates, screening caspase inhibitors, or as a positive control in caspase activity assays. We recommend using 1 unit/assay for analyzing caspase activity.

For a complete caspase assay protocol, please refer to BioVision's Caspase Fluorometric or Colorimetric Assay Kits.

**FOR RESEARCH USE ONLY! Not to be used in humans.**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- $\beta$ -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Quality Antibodies for Apoptosis and Signal Transduction Molecules