

## Caspase-Family Inhibitor, Boc-D-FMK

<b>CATALOG #:</b>	1120-20C
<b>STORAGE CONDITIONS:</b>	Store at -20°C.
<b>SHELF LIFE:</b>	1 year under proper storage conditions
<b>MOLECULAR WEIGHT:</b>	263.26
<b>SEQUENCE:</b>	Boc-Asp(OMe)-FMK (FMK, fluoromethyl ketone)
<b>PURITY:</b>	>99% by HPLC analysis.
<b>PACKAGE CONTENTS:</b>	20 µl (10 mM) of Boc-D-FMK in DMSO

### DESCRIPTION:

A synthetic peptide that irreversibly inhibits ICE-family protease/caspase activity and blocks apoptosis. The inhibitor is designed as a methyl ester to facilitate cell permeability. **(CAUTION:** If the intended use is on purified or recombinant enzymes, esterase should be added to generate free carboxyl groups).

### RECOMMENDED USAGE:

We recommend using 1000X dilutions for inhibiting caspase activity in Jurkat cell culture (e.g., Add 1 µl to 5 ml of culture medium). The optimal doses may vary for different cells and culture conditions.

**FOR RESEARCH PURPOSE ONLY! Not to be used in humans!**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- Cloning Insert Quick Screening Kit 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Monoclonal and Polyclonal Antibodies