

Caspase-4 Inhibitor, Z-LEVD-FMK

CATALOG #:	1144-1	1 mg
	1144-5	5 mg
STORAGE:	Store at -20°C.	
SHELF LIFE:	1 year under proper storage conditions	
MOLECULAR WEIGHT:	652.0	
SEQUENCE:	Z-Leu-Glu(OMe)-Val-Asp(OMe)-FMK	
PURITY:	One spot by TLC	
PACKAGE CONTENTS:	1 mg (or 5 mg) lyophilized powder	
RECONSTITUTION:	Reconstitute in DMSO (20 mM)	

DESCRIPTION:

A synthetic peptide that irreversibly inhibits activity of LEVD-dependent caspases (e.g., caspase-4). The inhibitor is designed as a methyl ester to facilitate cell permeability. The product can be used for both in vitro and in vivo studies. (**CAUTION:** If the intended use is on purified or recombinant enzymes, esterase should be added to generate free carboxyl groups).

RECOMMENDED USAGE:

We recommend using 2-10 μ M for inhibiting caspase activity in cell culture. The optimal doses may vary for different cells and culture conditions.

FOR RESEARCH PURPOSE ONLY! Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction Molecules