

96-well Plates (2)

CATALOG #: 1996-2

AMOUNT: 2 plates


INTRODUCTION:

These are polystyrene plates with clear C-bottom single-break strips provided as twelve 8-well strips in a support frame (12 x 8 matrix). The individual wells can be broken off separately (Single-break option). The C-bottom means a flat bottom with rounded corners. The flat bottom allows precise optics and the rounded corners allow pipetting without leaving a residue. The plates also offer a high-binding surface.

STORAGE CONDITIONS: Room temperature.

APPLICATIONS:

These plates are useful for diagnostic and immunological research applications and ideal for colorimetric assays using a 96-well format and standard plate readers.



FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

- 10K Spin Column
- 96-Well Plate 2 plate
- Avidin CovaLinked 96-well Plates
- Dounce Tissue Grinder
- Magnetic Separator (Re-usable)
- Protein G Coated 96-well Plates
- Metabolic Assay Kits
- Caspase Assay Kits

Features:

- (LxWxH): 127.8mm x 85.5 mm x 14.5 mm
- Strips mounted in frame
- C8 clear strips
- High binding surface