

FeTPPS

ALTERNATE NAME: 5,10,15,20-Tetrakis(4-sulfonatophenyl)porphyrinato Iron (III), Chloride

CATALOG #: 2065-25, 100

AMOUNT: 25 mg, 100 mg

STRUCTURE:



MOLECULAR FORMULA: $C_{44}H_{28}ClFeN_4O_{12}S_4$

MOLECULAR WEIGHT: 1024.27

CAS NUMBER: 90384-82-0

APPEARANCE: Dark brown solid

SOLUBILITY: H_2O

PURITY: $\geq 95\%$ by TLC

STORAGE: Store at $-20\text{ }^\circ\text{C}$. Protect from air and light

DESCRIPTION: A peroxyxynitrite decomposition catalyst. An iron porphyrin complex that catalytically isomerizes peroxyxynitrite to nitrate both in vivo and in vitro thus serving as a decomposition catalyst and a selective peroxyxynitrite scavenger. Does not complex with nitric oxide and exhibits little SOD mimetic activity.

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

RELATED PRODUCTS:

- Ascorbic Acid Assay Kit (**Cat. No. K661-100**)
- Ascorbic Acid Assay Kit II (FRASC) (**Cat. No. K671-100**)
- Caffeic Acid (**Cat. No. 2006-100, 500**)
- Catalase Assay kit (**Cat. No. K773-100**)
- Deferoxamine Mesylate (**1883-500, 1000**)
- (-)-Epigallocatechin gallate (**Cat. No. 1841-50**)
- Glutathione (GSH/GSSG Total) Assay Kit (**Cat. No. K264-100**)
- Glutathione Colorimetric Assay Kit (**Cat. No. K261-100**)
- Glutathione Fluorometric Assay Kit (**Cat. No. K251-100**)
- Glutathione Peroxidase Assay Kit (**Cat. No. K762-100**)
- Glutathione Reductase Assay Kit (**Cat. No. K761-100**)
- GST Colorimetric Assay Kit (**Cat. No. K263-100**)
- GST Fluorometric Assay Kit (**Cat. No. K263-100**)
- Hydrogen Peroxide Assay Kit (**Cat. No. K265-100**)
- Luteolin (**Cat. No. 2005-5, 25**)
- MnTBAP Chloride (**Cat. No. 2058-10, 50**)
- MnTMPyp Pentachloride (**Cat. No. 2056-10, 50**)
- Nitric Oxide Colorimetric Assay Kit (**Cat. No. K263-100**)
- Nitric Oxide Fluorometric Assay Kit (**Cat. No. K252-100**)
- Piperlongumine (**Cat. No. 1919-10, 50**)
- Piperlonguminine (**Cat. No. 1988-10, 50**)
- Resveratrol (**Cat. No. 1758-100, 500**)
- Rosmarinic Acid (**Cat. No. 1843-50, 250**)
- SOD Activity Assay Kit (**Cat. No. K335-100**)
- Thioredoxin Reductase Assay Kit (**Cat. No. K763-100**)
- Tiliroside (**Cat. No. 1879-1,5**)
- Total Antioxidant Capacity Assay Kit (**Cat. No. K274-100**)
- Uric Acid Assay Kit (**Cat. No. K608-100**)
- Xanthine Oxidase Assay Kit (**Cat. No. K710-100**)

FOR RESEARCH USE ONLY! Not to be used on humans.