

PRODUCT: **Tris-HCl Solution (1 M, pH 6.8)**

ALTERNATE NAME:	Tris hydrochloride
CATALOG #:	2106-100
AMOUNT:	100 ml
MOLECULAR WEIGHT:	157.60
CAS NUMBER:	1185-53-1
APPEARANCE:	Liquid
CONCENTRATION:	1 M
pH:	6.8
STORAGE:	Room temperature
DESCRIPTION:	A 1 M solution of Tris-HCl in DI water (pH 6.8)

RELATED PRODUCTS:

- Magnesium Chloride (1 M) **(2105)**
- Potassium Chloride Solution (1 M) **(2115)**
- 10% SDS in Water **(2102)**
- Sodium Chloride Solution (5 M) **(2101)**
- Sodium Chloride (0.9%), Sterile-Filtered **(9628)**
- Tris-HCl (pH 6.8) **(2106)**
- Tris-HCl (pH 8.8) **(2107)**
- Triton X-100 **(2104)**
- Tween 20 (10%) **(2109)**
- Wash Buffer **(1210)**

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

USAGE: **FOR RESEARCH USE ONLY! Not to be used in humans**