


PRODUCT: AdipoRon hydrochloride**ALTERNATE NAME:** 2-(4-Benzylphenoxy)-N-(1-benzylpiperidin-4-yl)acetamide hydrochloride**CATALOG #:** 2572-5, 25**AMOUNT:** 5 mg, 25 mg**STRUCTURE:****MOLECULAR FORMULA:** C₂₇H₂₈N₂O₃·HCl**MOLECULAR WEIGHT:** 464.98**CAS No.** 1781835-20-8**APPEARANCE:** White solid**SOLUBILITY:** DMSO (>10 mg/ml) or EtOH (>5 mg/ml) or H₂O (1 mg/ml)**PURITY:** ≥98% by HPLC and NMR**STORAGE:** Store at -20°C. Protect from air and moisture**DESCRIPTION:** A water-soluble derivative of the potent and selective adiponectin receptor agonist AdipoRon (Cat. No. 2565).**REFERENCES:** Okada-Iwabu, M., *et al.* (2013). *Nature* **503**, 493-499.**RELATED PRODUCTS:**

Adiponectin, human recombinant (Cat. No. 4901-25, -100, -1000)
Adiponectin, murine recombinant (Cat. No. 4902-25, -100, -1000)
Adiponectin, rat recombinant (Cat. No. 4903-25, -100, -1000)
gAcrp30/Adiponectin, human recombinant (Cat. No. 4591-25, -100, -1000)
gAcrp30/Adiponectin, murine recombinant (Cat. No. 4592-25, -100, -1000)
gAcrp30/Adiponectin, rat recombinant (Cat. No. 4593-25, -100, -1000)
Adiponectin Antibody (Cat. No. 5901-50, -250)
Adiponectin Antibody (Cat. No. 5902-50, -250)
Adiponectin Antibody (Cat. No. 5903-50, -250)
Adiponectin Antibody (Clone MAD1-1147) (Cat. No. 5900-50)
Adiponectin (human) Elisa Assay Kit (Cat. No. K4901-100)
Adiponectin (mouse) Elisa Assay Kit (Cat. No. K4902-100)
Adiponectin (rat) Elisa Assay Kit (Cat. No. K4903-100)
AdipoRon (Cat. No. 2565-5, 25)

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.**USAGE:** **FOR RESEARCH USE ONLY! Not to be used in humans.**