

Bax Polyclonal Antibody

CATALOG NO: 3032-30T 30 µg (Trial size)
3032-100 100 µg

HOST: Rabbit

IMMUNOGEN: Synthetic peptide corresponding to amino acids near the N-terminus of human Bax (ID# BV-126).

FORMULATION:

0.5 mg/ml affinity purified rabbit anti-Bax polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 30% glycerol, 0.5% BSA, 0.01% thimerosal.

SPECIES REACTIVITY: Human, Mouse, Rat, Monkey

STORAGE CONDITIONS:

Store at -20° C. For long term storage, aliquot and freeze at -70° C. Avoid repeated freeze/thaw cycles.

DESCRIPTION:


The Bcl-2 (B-cell leukemia/lymphoma 2) family proteins function to either inhibit (Bcl-2 and Bcl-X_L) or promote (Bax and Bak) apoptosis. Bax contains the conserved Bcl-2 homology 1 (BH1) and 2 (BH2) domains that allow for its homodimerization or heterodimerization with Bcl-2. When Bcl-2 is in excess, apoptosis is inhibited. However, if Bax levels increase in response to a death signal, the cell is pushed toward death.

SPECIFICITY:

The antibody detects a 20 kDa protein, corresponding to the apparent molecular weight of Bax on SDS-PAGE immunoblots. Rat kidney tissue lysate can be used as a positive control.

APPLICATION AND USAGE:

The antibody can be used for Western blot analysis (0.5-4 µg/ml). Based on researchers feedback, it can also be used in immunoprecipitation (10-20 µg/ml). However, the optimal conditions should be determined individually.



FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Cholesterol & HDL/LDL Quantification Kits

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies