

## Anti-Rat TNF- $\alpha$ Polyclonal Antibody

<b>CATALOG NO:</b>	3054-30T	30 $\mu$ g (Trial size)
	3054-100	100 $\mu$ g
<b>HOST:</b>	Rabbit	
<b>IMMUNOGEN:</b>	Recombinant rat TNF- $\alpha$ (Internal ID# BV-A23-1)	
<b>SPECIES:</b>	Rat	

### FORMULATION:

0.5 mg/ml affinity purified polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 30% glycerol, 0.5% BSA, 0.01% thimerosal.

### DESCRIPTION:


Tumor Necrosis Factor-alpha (TNF- $\alpha$ ) is a potent lymphoid factor that exerts cytotoxic effects on a wide range of tumor cells and certain other target cells. Rat TNF- $\alpha$  is a 17 kDa protein containing 157 amino acid residues.

### STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

### APPLICATION AND USAGE:

The antibody can be used in Western blot analysis (0.5-4  $\mu$ g/ml). Recombinant rat TNF- $\alpha$  (Cat.# 1052-10) and rat brain tissue lysate can be used as positive controls. Per researcher's feedback, the antibody can also be used in immunoprecipitation and Neutralization (2-10  $\mu$ g/ml). However, the optimal conditions should be determined individually.



### Cat. #:3054-100

Western blot analysis with recombinant Rat TNF- $\alpha$ .

Lane 1: 10 ng

Lane 2: 50 ng

Lane 3: 250 ng

Rat TNF- $\alpha$

**FOR RESEARCH USE ONLY! Not to be used on humans.**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- $\beta$ -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Monoclonal and Polyclonal Antibodies