

ASK1/MAPKKK5 Polyclonal Antibody

CATALOG #: 3128-100

LOT #: _____

HOST: Rabbit

IMMUNOGEN: Synthetic peptide corresponding to residue surrounding amino acid 1359 of human ASK1.

FORMULATION:

100 µg (0.2 mg/ml) immunoaffinity purified, rabbit anti-ASK1 polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA, 0.02% sodium azide.

SPECIES REACTIVITY: Human

STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and refreeze at -70°C. Avoid repeated freeze/thaw cycles.

DESCRIPTION:

ASK1 (Apoptosis Signal-regulating Kinase 1) activates two different subgroups of MAPKK, MKK4 and MKK6, which in turn activate c-Jun N-terminal Kinase (JNK) and p38 MAP Kinase. ASK1 is activated by TNFR and Fas through the interaction with members of the TRAF Family and the Fas-associated protein Daxx. Overexpression of ASK1 induces apoptosis. However, a catalytically inactive form of ASK1 has been shown to inhibit TNF-α induced apoptosis.

SPECIFICITY:

The immunoaffinity purified antibody detects the 155 kDa ASK1 protein in human SW1353 cell lysate by immunoblotting.

APPLICATION AND USAGE:

The antibody can be used for Western blot analysis (1-4 µg/ml). However, the optimal conditions should be determined individually.

FOR RESEARCH USE ONLY! Not to be used in human subjects!

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies