

## NIK Polyclonal Antibody

<b>CATALOG NO:</b>	3193R-30T	30 µg (Trial size)
	3193R-100	100 µg
<b>HOST:</b>	Rabbit	
<b>IMMUNOGEN:</b>	Synthetic peptide surrounding amino acid 927 of human NIK (Internal ID# BV-F49G)	
<b>SPECIES REACTIVITY:</b>	Human, Mouse, Rat, Horse, dog	

### FORMULATION:

0.5 mg/ml affinity purified rabbit anti-NIK polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 30% glycerol, 0.5% BSA, 5 mM EDTA and 0.03% proclin.

### STORAGE CONDITIONS:


Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/defrost cycles.

### BACKGROUND DESCRIPTION:

NIK (NF-κB-inducing kinase) is a member of the MAP kinase kinase kinase family that binds TRAF2 and stimulates NF-κB activity. NIK was initially isolated from a human B cell cDNA library and contains 795 amino acids with an apparent molecular weight of slightly more than 97 kDa on SDS gel. NIK is a serine/threonine kinase and its kinase activity contributes to IκB phosphorylation. The carboxyl terminal segment of NIK binds TRAF2. A mutant NIK with intact carboxyl terminus but without the two lysine residues at its catalytic domain serves as a dominant-negative inhibitor for NF-κB activation. NIK also interacts with TRAF6 and mediates IL-1-induced NF-κB activation.

### APPLICATION AND USAGE:

The antibody can be used for Western blotting (0.5-4 µg/ml). However, the optimal conditions should be determined individually. Blocking peptide (Cat.# 3191RRBP-50) is available separately.



**Cat. #: 3193R-100**  
Western blot analysis of NIK using Jurkat cell lysate (Lane 1 & 2)

**FOR RESEARCH USE ONLY! Not to be used on humans.**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Monoclonal and Polyclonal Antibodies