

ACADVL Polyclonal Antibody

CATALOG NO:	3240-30T	30 µg (Trial size)
	3240-100	100 µg
IMMUNOGEN:	Synthetic peptide surrounding amino acid 441 of rat ACADVL	
INTERNAL ID:	BV-K81	
HOST:	Rabbit	
SPECIES REACTIVITY:	human, mouse, rat	

FORMULATION:

0.5 mg/ml affinity purified rabbit anti-ACADVL polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 30% glycerol, 0.5% BSA, 5 mM EDTA 0.5% BSA, 5 mM EDTA and 0.03% proclin.

STORAGE CONDITIONS:


Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

DESCRIPTION:

Acyl-CoA dehydrogenase, very long chain (ACADVL) is a member of the acyl-CoA dehydrogenases family. ACADVL is an inner mitochondrial membrane protein that catalyzes the first step of the mitochondrial fatty acid beta-oxidation pathway. A deficiency of ACADVL reduces myocardial fatty acid beta-oxidation and is associated with cardiomyopathy.

APPLICATION AND USAGE:

The antibody can be used for Western blot analysis (1-4 µg/ml). However, the optimal conditions should be determined individually. Blocking peptide (Cat.# 3240BP-50) is available separately.



ACADVL

Cat. #: 3240-100 ACADVL Antibody.
Western blot analysis of ACADVL with rat kidney tissue lysate.

FOR RESEARCH USE ONLY! Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies