

Phospho-AFX (Ser193) Polyclonal Antibody

CATALOG #: 3270-100

AMOUNT: 100 µg

LOT #: _____

FORMULATION:

100 µg (0.5 mg/ml) peptide affinity purified polyclonal rabbit anti-Phospho-AFX antibody in phosphate-buffered saline (PBS) containing 50% glycerol, 1% BSA, and 0.02% sodium azide

SPECIES REACTIVITY: Human

STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

DESCRIPTION:

The forkhead family of transcription factors is involved in tumorigenesis in rhabdomyosarcoma and acute leukemias. Within the family, three members (FKHR, AFX, AND FKHL1) have sequence similarity to the nematode homolog Daf 16, recently shown to mediate signaling via a pathway involving IGFR1, P13-K and Akt. There are three potential Akt phosphorylation sites on the FKHR proteins. One site, FKHR Ser256 is found in the Forkhead domain and the other two are outside this domain. Recent studies show that phosphorylation of FKHR family members by Akt promotes cell survival. Phosphorylation regulates FKHR nuclear translocation and target gene transcription. Insulin also stimulates FKHR phosphorylation at Ser253.

SPECIFICITY:

Phospho-AFX (Ser193) antibody detects only the phosphorylated Ser193 of AFX and does not react with nonphosphorylated AFX.

APPLICATION AND USAGE:

The antibody can be used for Western blotting (1-2 µg/ml). However, the optimal concentrations should be determined individually.

FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies