

## Heme Oxygenase-2 Polyclonal Antibody

<b>CATALOG NO:</b>	3392-30T      30 µg (Trial size) 3392-100      100 µg
<b>HOST:</b>	Rabbit
<b>IMMUNOGEN:</b>	Synthetic peptide surrounding amino acid 34 of human HO-2 (Internal ID# BV-A2)
<b>SPECIES REACTIVITY:</b>	Human, Mouse, Rat, Hamster, Monkey

### FORMULATION:

0.5 mg/ml affinity purified rabbit polyclonal antibody in phosphate-buffered saline (PBS) containing 30% glycerol, 0.5% BSA, and 0.01% thimerosal.

### STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

### BACKGROUND:


Heme oxygenase (HO) is a microsomal enzyme that oxidatively cleaves heme, a pro-oxidant, into carbon monoxide and biliverdin. Two heme oxygenase isozymes, termed HO-1 and HO-2, have been identified so far. At the amino acid level, HO-1 and HO-2 are 42% homologous. However, each isozyme is evolutionary highly conserved. HO-1 is the inducible form of the enzyme and is therefore also referred to as stress protein Hsp32. Expression of HO-1 can be induced by a variety of stimuli such as heme, metals, hormones, UV radiation and sulfhydryl depleting agents. In contrast, HO-2 is the constitutive form of the enzyme. It is not inducible by exogenous stimuli and HO-2 is the most prevalent form in most tissues, except the spleen where HO-1 levels are predominant under normal conditions.

### SPECIFICITY:

The antibody recognizes ~36 kDa Heme-Oxygenase-2 of human, mouse, rat, hamster, and monkey origins. Blocking peptide (Cat.# 3392BP-50) is available separately.

### APPLICATION:

The antibody has been used in Western blotting (0.5-4 µg/ml). Per researcher's feedback, it can also be used in immunoprecipitation (20 µg/ml) and immunohistochemistry (20 µg/ml, frozen section tested). However, the optimal concentrations should be determined individually.



### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

**FOR RESEARCH USE ONLY! Not to be used on humans.**