

HNF-4 Polyclonal Antibody

CATALOG NO:	3688-30T	30 µg (Trial size)
	3688-100	100 µg
HOST:	Rabbit	
IMMUNOGEN:	Synthetic peptide surrounding amino acid 454 of rat HNF (ID# BV-242).	
SPECIES REACTIVITY:	Mouse, Rat	

FORMULATION:

0.5 mg/ml affinity purified rabbit polyclonal antibody in phosphate-buffered saline (PBS) containing 30% glycerol, 0.5% BSA, and 0.01% thimerosal.

STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw.

DESCRIPTION:


HNF-4 (Hepatocyte nuclear factor 4), a member of the steroid/thyroid nuclear receptor superfamily, is a transcriptional factor which expresses in the liver, intestine, kidney, and pancreatic cells. It contains several functional domains: a ligand-independent activation domain (AF1), a zinc finger DNA binding domain, and a ligand-dependent activation domain (AF2). HNF-4 binds to a specific DNA element as a homodimer and regulates the expression of many genes, involved in glucose, fatty acid, and cholesterol metabolisms.

SPECIFICITY:

The antibody recognizes HNF-4 from samples of mouse and rat origins. Reactivity to other species has not been determined.

APPLICATION:

The antibody can be used for Western blotting (0.5-4 µg/ml). However, the optimal concentrations should be determined individually. Mouse small intestine and rat kidney (Cat.# 2702-100) tissue lysates can be used as positive controls.



FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies