

## Recombinant Human BAFF

<b>CATALOG #:</b>	4451-20	20 µg
	4451-100	100 µg
	4451-1000	1 mg
<b>LOT #:</b>	_____	
<b>SOURCE:</b>	<i>E. coli</i>	
<b>PURITY:</b>	>95% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of BAFF.	
<b>MOL. WEIGHT:</b>	19.335 kDa	
<b>FORM:</b>	Lyophilized from 0.3X PBS, pH 7.5	

**RECONSTITUTION:**

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. This solution can then be diluted into other aqueous buffers or stored at 4°C for 1 week or – 20°C for future use.

**STORAGE CONDITIONS:**

The lyophilized BAFF is best-stored desiccated below 0°C. Reconstituted human BAFF should be stored in working aliquots at -20°C.

**DESCRIPTION:**

Human BAFF (for B cell activating factor belonging to the TNF family) is a newly discovered novel ligand of the TNF family. BAFF plays an important role as co-stimulator of B cell proliferation and function. Recombinant human BAFF is a soluble protein containing 153 amino acid residues + an N-terminal His-tag (MW 19.335 kDa).

**BIOLOGICAL ACTIVITY:**

Determined by the dose-dependent stimulation of IL-8 production by human PBMC. The ED<sub>50</sub> for this effect is less than 10 ng/ml. Please note: Results may vary with PBMC donors.

**RELATED PRODUCTS:**

## Apoptosis Detection Kits &amp; Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

## Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

## Cell Proliferation &amp; Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

## Cell Damage &amp; Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

## Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

## Adipocyte &amp; Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

## Molecular Biology &amp; Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

## Growth Factors and Cytokines

## Monoclonal and Polyclonal Antibodies

**FOR RESEARCH USE ONLY! Not to be used in humans.**