

Recombinant Murine PDGF-AA

CATALOG #:	4483-10	10 µg
	4483-1000	1 mg
LOT #:	_____	
SOURCE:	<i>E. coli</i>	
PURITY:	> 97% by SDS-PAGE and RP-HPLC analysis Endotoxin level is < 0.01 ng per µg of rm-PDGF-AA.	
MOL. WEIGHT:	28.9 kDa	
FORM:	Lyophilized with no additives.	

RECONSTITUTION:

Centrifuge vial before opening. When reconstituting the product, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution. It is recommended to reconstitute the lyophilized product with sterile water to a concentration of 0.1-0.5 mg/ml, which can be further diluted into other aqueous solutions.

STORAGE CONDITIONS:

The lyophilized rm-PDGF-AA is very stable at -20 °C. Upon reconstitution the product should be aliquoted and stored at -20 °C. It is recommended to add a carrier protein (0.1% HSA or BSA) for long term storage.

DESCRIPTION:

Recombinant Mouse Platelet Derived Growth Factor-AA is a disulfide linked dimer comprised of 2 polypeptide chains of 126 amino acids each, with a dimeric molecular weight of 28.9 kDa. Belonging to a family of PDGF's including PDGF-AA, PDGF-BB and PDGF-AB, they play an important role as a mitogen for a number of cell types. Recombinant Mouse PDGF-AA produced in E.Coli is a non-glycosylated homodimer.

BIOLOGICAL ACTIVITY:

The activity is determined by the dose-dependent proliferation of Balb/c 3T3 cells and is typically 8 -10 ng/ml.

AA Sequence:

MSIEEAVPAV CKTRTVIYEI PRSQVDPTSA NFLIWPPCVE VKRCTGCCNT SSVKQCPSRV
HHRSVKVAKV EYVRKKPKLK EVQVRLEEHL ECACATSNLN PDHREEETGR RRESGKNRKR
KRLKPT

FOR RESEARCH USE ONLY! Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies