

Recombinant Human IL-22

CATALOG #:	4601-10	10 µg
	4601-50	50 µg
	4601-1000	1 mg
SOURCE:	<i>E. coli</i>	
PURITY:	>97% by SDS-PAGE and HPLC Endotoxin level is <0.1 ng per µg of IL-22.	
MOL. WEIGHT:	16.8 kDa	
FORM:	Lyophilized from 1 mM Sodium Citrate, pH 4.0.	

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. It is recommended that further dilutions be made into buffer containing carrier protein or medium containing serum. This solution can then be stored at 4°C for 1 week or -20° for future use.

STORAGE CONDITIONS:

The lyophilized IL-22 is best-stored desiccated below 0°C. Reconstituted IL-22 should be stored at working aliquots at -20°C.

DESCRIPTION:

Human Interleukin-22 (IL-22), also known as T cell-derived inducible factor (ILTIF), is a cytokine that is structurally related to IL-10. Studies have shown that IL-22 induces acute-phase reactants *in vitro* and *in vivo*, suggesting that the cytokine contributes to the regulation of inflammatory responses. Recombinant human IL-22 is a 16.8 kDa protein consisting of 147 amino acids residues.

ACTIVITY:

Recombinant IL-22 activates STAT following receptor-ligand interaction.

FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies