

## Recombinant Murine VEGF-C

<b>CATALOG #:</b>	4634-10	10 µg
	4634-50	50 µg
	4634-1000	1 mg
<b>LOT #:</b>	_____	
<b>SOURCE:</b>	<i>E. coli</i>	
<b>PURITY:</b>	>90% by SDS-PAGE analyses Endotoxin level is <0.1 ng per µg of VEGF-C	
<b>MOL. WEIGHT:</b>	15 kDa	
<b>FORMULATION:</b>	Lyophilized without additives	

**RECONSTITUTION:**

Centrifuge the vial prior to opening. Reconstitute in 0.1% acetic acid to a concentration of 0.1-1.0 mg/ml. The solution should be stored at -70°C.

**STORAGE CONDITIONS:**

The lyophilized VEGF-C is best-stored at -20°C. Reconstituted VEGF-C should be stored in working aliquots at -70°C.

**DESCRIPTION:**

VEGF-C (Vascular Endothelial Growth Factor-C) is a 125 amino acid recombinant protein that lacks the N-terminal and C-terminal extensions and contains only the middle third VEGF homology domain.

**ACTIVITY:**

Testing in progress.

**RELATED PRODUCTS:**

## Apoptosis Detection Kits &amp; Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

## Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

## Cell Proliferation &amp; Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

## Cell Damage &amp; Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

## Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

## Adipocyte &amp; Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

## Molecular Biology &amp; Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

## Growth Factors and Cytokines

## Monoclonal and Polyclonal Antibodies

**FOR RESEARCH USE ONLY! Not to be used on humans.**