

## Recombinant Human VEGF-B

<b>CATALOG #:</b>	4642-20	20 µg
	4642-100	100 µg
	4642-1000	1 mg
<b>LOT #:</b>	_____	
<b>SOURCE:</b>	<i>E. coli</i>	
<b>PURITY:</b>	>98% by SDS-PAGE analyses Endotoxin level is <0.1 ng per µg of VEGF-B	
<b>MOL. WEIGHT:</b>	19.63 kDa	
<b>FORM:</b>	Lyophilized without additives	

### RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in 0.1% acetic acid to a concentration of 0.1-1.0 mg/ml. The solution should be stored at -70°C.

### STORAGE CONDITIONS:

The lyophilized VEGF-C is best-stored at -20°C. Reconstituted VEGF-C should be stored in working aliquots at -70°C.

### DESCRIPTION:

VEGF-B, a member of the VEGF family, is a potent growth and angiogenic cytokine. It promotes DNA synthesis in endothelial cells, regulates angiogenesis and vascular permeability, and inhibits apoptosis in certain smooth muscle cells and neurons. VEGF-B is expressed in all tissues except the liver. It forms cell surfaced-associated disulfide linked homodimers and can form heterodimers with VEGF-A. There are two known isoforms, formed by alternative splicing, which have been designated VEGF-B167 and VEGF-B186. Both forms have identical amino-terminal sequences encoding a "cysteine knot" like structural motif, but differ in their carboxyl-terminal domains. Both VEGF-B isoforms signal only through the VEGFR1 receptor. Recombinant human VEGF-B is a 38.0 kDa disulfide-linked homodimeric protein consisting of two 167 amino acid polypeptide chains.

### BIOLOGICAL ACTIVITY:

Per researcher's feedback, the protein is biologically active as determined by the dose-dependent stimulation of the proliferation of human umbilical vein endothelial cells (HUVEC) in the presence of human VEGF<sub>165</sub>. The expected ED<sub>50</sub> for this effect is 0.5-5.0 µg/ml.

**FOR RESEARCH USE ONLY! Not to be used on humans.**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Monoclonal and Polyclonal Antibodies