

Recombinant Human Gremlin-2

CATALOG #:	4657-20	20 µg
	4657-100	100 µg
	4657-1000	1 mg
LOT #:	_____	
SOURCE:	<i>E. coli</i>	
PURITY:	≥95% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of Gremlin-2	
MOL. WEIGHT:	17.2 kDa	
FORMULATION:	Lyophilized with no additives	

RECONSTITUTION:


Centrifuge the vial prior to opening. Reconstitute to 1 mg/ml with sterilized dH₂O. Pipette to dissolve the protein pellet completely.

STORAGE CONDITIONS:

The lyophilized protein is best-stored desiccated below 0°C. Reconstituted protein should be stored at working aliquots at -20°C. Avoid freeze/thaw cycles.

DESCRIPTION:

Gremlin (also known as Increased in High Glucose Protein 2, IHG-2, Down-regulated in Mos-transformed cells protein, Drm) functions as a bone morphogenetic protein (BMP) antagonist. It acts by binding to, and forming heterodimers with BMP-2, BMP-4, BMP-7, thus preventing them from interacting with their cell surface receptors. Recombinant human Gremlin-2 is a N-terminus His-tagged protein that was expressed from *E. coli*.



Cat. #: 4657-20

Purity of recombinant human Gremlin-2 (2 µg) was analyzed by SDS-PAGE electrophoresis.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

FOR RESEARCH USE ONLY! Not to be used on humans.