

## Recombinant Mouse PLGF

<b>CATALOG #:</b>	4743-10	10 µg
	4743-50	50 µg
	4743-1000	1 mg
<b>LOT #:</b>	_____	
<b>SOURCE:</b>	<i>E. coli</i>	
<b>PURITY:</b>	>95% by SDS-PAGE and HPLC analysis	
<b>ENDOTOXIN:</b>	Endotoxin level is <0.1 ng per µg of PLGF.	
<b>MOL. WEIGHT:</b>	32.6 kDa (2 x16.3 kDa disulfide-linked homodimer)	
<b>FORMULATION:</b>	Sterile filtered and lyophilized without additives	

**RECONSTITUTION:**

Centrifuge the vial prior to opening. It is recommended to reconstitute the protein in sterile 20 mM acetic acid to a concentration of 0.1-1.0 µg/µl. The solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

**STORAGE CONDITIONS:**

The lyophilized protein is best-stored desiccated below 0°C. Reconstituted protein should be stored at working aliquots at -20°C. For long term storage it is recommended to add a carrier protein (0.1% HAS or BSA). Prevent freeze/thaw cycles.

**DESCRIPTION:**

PLGF (also called Placenta Growth Factor, PIGF, PGF) is an angiogenic factor that belongs to the cysteine-knot superfamily of growth factors. PLGF is expressed in placental tissues, colon and mammary carcinomas. It signals through the VEGFR-1/FLT1 receptor and stimulates endothelial cell proliferation and migration. Recombinant mouse PLGF is a disulfide-linked homodimer containing two 142 amino acid polypeptide chains.

**Biological Activity:**

Determined by its ability to chemoattract human monocytes of using a concentration range of 0.5-10 ng/ml.

**FOR RESEARCH USE ONLY! Not to be used in humans.**

**RELATED PRODUCTS:**

## Apoptosis Detection Kits &amp; Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

## Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

## Cell Proliferation &amp; Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

## Cell Damage &amp; Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

## Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

## Adipocyte &amp; Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

## Molecular Biology &amp; Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

## Growth Factors and Cytokines

## Monoclonal and Polyclonal Antibodies