

Recombinant Human Erythropoietin- α (EPO- α)

CATALOG #: 4764-500 500IU

SOURCE: Chinese Hamster Ovary Cells (CHO)

PURITY: \geq 98% as determined by SDS-PAGE and HPLC

ENDOTOXIN: Less than 0.1 ng/ μ g (IEU/ μ g) of rhEPO

PHYSICAL APPEARANCE: Lyophilized powder

FORMULATION: Each mg of lyophilized protein contains 0.58 mg sodium citrate, 0.58 mg sodium chloride and 0.006 mg citric acid.

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1 to 1.0 mg/ml. This solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

STORAGE CONDITIONS:

The lyophilized protein is best-stored desiccated below 0°C. Reconstituted human EPO should be stored at working aliquots at -20°C. Avoid freeze-thaw cycles.

DESCRIPTION:

Recombinant Human Erythropoietin-alpha produced in Chinese hamster ovary (CHO) cells by recombinant DNA technology is a single, polypeptide chain containing 165 amino acids and having a predicted molecular mass of 21 kDa and apparent glycosylated molecular mass of 30.4 kDa. EPO-alpha is purified by proprietary chromatographic techniques.

ACTIVITY:

BioVision's Erythropoietin is fully biologically active when compared to standard. The Specific Activity measured by Normocyth-aemic mice was found to be 120,000IU/mg.

USAGE: For Research Use Only. Not to be used in humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Additional Apoptosis Kits & Reagents

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione Fluorometric & Colorimetric Assay Kits
- Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- Camp & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit

Monoclonal & Polyclonal Antibodies (many)

Growth Factors & Cytokines (many)