

## Recombinant Human TFF1

**CATALOG #:** 4893-20 20 µg  
4893-100 100 µg  
4893-1000 1 mg

**LOT #:** \_\_\_\_\_

**SOURCE:** *E. coli*

**PURITY:** >98% by SDS-PAGE and HPLC analyses  
Endotoxin level is <0.1 ng per µg of human TFF1.

**MOL. WEIGHT:** 8.2 kDa

**FORM:** Sterile filtered and lyophilized powder

### RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. The solution can then be diluted into other aqueous buffers and store at 4°C for 1 week or –20°C for future use.

### STORAGE CONDITIONS:


The lyophilized protein is best-stored desiccated below 0°C. Reconstituted TFF1 should be stored at working aliquots at -20°C.

### DESCRIPTION:

The Trefoil Factor peptides (TFF1, TFF2 and TFF3) are expressed in the gastrointestinal tract, and appear to play an important role in intestinal mucosal defense and repair. Recombinant human TFF1 is a 8.2 kDa polypeptide containing 60 amino acid residues.

### BIOLOGICAL ACTIVITY:

Determined by its ability to chemoattract human MCF-7 cells using a concentration 5-10 µg/ml.



### Cat. # 4893-20

2 µg of rh-TFF-1 was subjected to SDS-PAGE and Coomassie blue staining. The molecular weight is ~8.2 kDa.

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction Molecules

**FOR RESEARCH USE ONLY! Not to be used on humans.**